

**Kiato** <sup>+</sup> PLUS


Sterile podiatric  
Gouge blades




**DIAMANT**




*The smart choice for sharp blades*

- Swiss carbon steel 
- Sterilized by Gamma Radiation to a dosage of 2.5 Mrads.
- Greater strength due to perfectly formed channel.
- Durable sharpness with smooth cutting.



Size	1	2	3	4	5	6	8	10	12	15	18
Ref. No.	12101	12102	12103	12104	12105	12106	12108	12110	12112	12115	12118

Produced on SWISS CNC Machine with vision control for better edge geometry



Individually foil packed in box of 50 blades & 100 blades. 10 boxes per case. 5000 pcs per master shipper carton.



MDSS GmbH  
Schiffgraben 41  
30175 Hannover  
Germany



Aditya Dispomed Products Pvt. Ltd.  
Plot No. 19 Sector 6 IMT Manesar  
Gurgaon Haryana 122050 India

Tel. +91 124 47 64 900  
Fax +91 124 47 64 911

info@adityadispmomed.com  
www.adityadispmomed.com

aditya  
DISPOMED

A subsidiary of Kehr Surgical Pvt. Ltd.

All products are CE marked as per the European Medical Device Directive 93/42/EEC.