

Fine Blades

STAINLESS STEEL SURGICAL BLADES




- Stainless steel blades made of premium Swedish Steel.
- Individually packed in easy peel foil Pack.
- Fits all makes of collet handles.
- Sterilized by Gamma Radiation to a minimum dosage of 2.5 Mrads.
- Blades confirm to Medical Device Directive 93/42/EEC.

The Fine Range - accuracy you demand and quality as standard

Complimenting the traditional blade shapes and developed potentially for more complex and challenging surgical procedures which may be encountered in such disciplines as reconstruction, ophthalmics, cardiology and even hair restoration, the surgeons may sometimes prefer our Fine Range.

The Fine Range has small, slim, precise blades in stainless steel for delicate operations.

Special handles, with independently screwed top sections to facilitate



All products are CE marked as per the European Medical Device Directive 93/42/EEC.

EC REP MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany

CH REP MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland



EDITION 011/1 June 23

aditya
DISPOMED

A subsidiary of Kehr Surgical Pvt. Ltd.



Aditya Dispomed Products Pvt. Ltd.
Plot No. 19 Sector 6 IMT Manesar
Gurgaon Haryana 122050 India

Tel. +91 124 22 91 406 / 7
Fax +91 124 22 91 206

info@adityadispomed.com
www.adityadispomed.com

KIATO⁺ PLUS

FINE STAINLESS STEEL SURGICAL BLADES



BB361R

A chisel shaped blade used for nail treatments in Podiatry and Chiropody.



BB362R

A chisel shaped blade used for nail treatments in Podiatry and Chiropody.



BB363R

Precisely ground on four sides like an arrow head, the BB363R can be used in a number of specialist disciplines including Plastic & Reconstructive Surgery for procedures such as Otoplasty, Rhytidectomy and Rhinoplasty.




BB364R

The BB364R can be used in a number of specialist disciplines including Plastic & Reconstructive Surgery, Hip and Knee Arthroscopy, Tenotomy and minimal invasive foot and ankle procedures.



BB365R

The BB365R can be used in a number of specialist disciplines including Plastic & Reconstructive Surgery, Orthopaedics, Cardio-Vascular Surgery and Dental Periodontics.



BB366R

The BB366R can be used in a number of specialist disciplines including Plastic & Reconstructive Surgery, Orthopaedics, Cardio-Vascular Surgery and Dental Periodontics.



BB367R

The BB367R can be used in a number of specialist disciplines including Plastic & Reconstructive Surgery, Ophthalmology, Cardio-Vascular Surgery and Dental Periodontics



BB368R

The distinctive hooked shaped BB368R can be used in a number of specialist disciplines including cosmetic procedures such as a "suction lipectomy", E.N.T. Surgery for "Tonsillectomy" and also in Dental Periodontics.



BB369R

Precisely ground on four sides like an arrow head, the BB369R can be used in a number of specialist disciplines including Plastic & Reconstructive Surgery, Cardiology, Ophthalmology, Endoscopy, Laparoscopy, Pyloromyotomy for infantile pyloric stenosis and Dental Implantology techniques.