


Kiato



The smart choice for sharp blades



Better sharpness. Better grip. Better feel. Better performance.

- Swedish Stainless Steel Blade.
- ABS plastic handles ensure greater strength.
- Graduated scale for convenience.
- Blade size printed on handle



Also offered in blade size 10A, 10R, 11P, 12D, 13, 18, 19, 25 & 26

All products are CE marked as per the European Medical Device Directive 93/42/EEC.



MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany



MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

10

11

12

15

15C

15


20

21

22

23

24



Mini Scalpel – Also known as thumb scalpel. Fitted with a stainless steel blade and safety plastic guard. Sold in bulk. An ideal product for kits.



Packed in easy peel pouches. 10 scalpels per unit box.
Sterilized by Ethylene Oxide



Boxes are colour coded for easy identification.
10 boxes per case. (Mini scalpels sold in bulk only)

⚠ Disposable scalpels are sharp. Handle with care. Do not touch the cutting edge. Dispose in an approved manner.



Aditya Dispomed Products Pvt. Ltd.
Plot No. 19 Sector 6 IMT Manesar
Gurgaon Haryana I22050 India
Tel. +91 124 47 64 900
Fax +91 124 47 64 911
info@adityadisposed.com
www.adityadisposed.com

A subsidiary of Kehr Surgical Pvt. Ltd.