


Kiato



The smart choice for sharp blades

Kiato

Sterile Disposable Scalpel

"Consistent sharpness. Secure grip. Reliable performance."

- Swedish Stainless Steel Blade.
- ABS plastic handles ensure greater strength.
- Graduated scale for convenience.
- Protective Cap



Also offered in blade size 10A, 10R, 11P, 12D, 13, 18, 19, 25 & 26

CE marked medical device.
Classified as a medical device under Regulation (EU) 2017/745 (MDR)



MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany



MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

- Single-use only
- Do not re-sterilize
- Do not use if package is damaged
- Dispose of used scalpel in an approved sharps container



Packed in easy peel pouches. 10 scalpels per unit box.

Sterilized by Ethylene Oxide



Size No colour coded labels affixed to the box for easy Identification. 10 scalpels per box. 10 boxes per inner case.

⚠ "Disposable scalpels are sharp. Handle with care. Avoid contact with the cutting edge. Dispose of after use in accordance with local sharps disposal regulations."

aditya
DISPOMED



Aditya Dispomed Products Pvt. Ltd.
Plot No. 19 Sector 6 IMT Manesar
Gurgaon Haryana 122050 India
Tel. +91 124 47 64 904

sales@adityadispomed.com
www.adityadispomed.com

The smart choice for sharp blades