



Scalpel Handle

Material

Stainless Steel X6 Cr 17 (1.4016)
Grade e as per ISO 7153-1 : 2016

Fitment Dimension

The slot conforms
to EN ISO 27740 (ISO 7740)

Carbon (C)	Silicon (Si)	Chromium (Cr)	Manganese (Mn)	Phosphorus (P)	Sulphur (S)
≤ 0.08%	≤ 1.00%	16.00%- 18.00%	≤ 1.00%	≤ 0.04%	≤ 0.015%

The **KIATO** metal BP Handle is made using German raw material, cold forged and precision machined for perfectly reliable blade fitment with no wobble.

Fitment No.	REF	Suitable for blade Sizes
3	10003	10, 10A, 10R, 11, 11P, 12, 12D, 13, 14, 15, 15T, 15c, 15SM
4	10004	18, 19, 20, 21, 22, 23, 24, 25, 26 & 36

Packaging

- Individually packed in a rigid blister
- 10 handles to a dispenser box

EAN Code : 08902895100037



MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany



MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland



Aditya Dispomed Products Pvt. Ltd.
Plot No. 19 Sector 6 IMT Manesar
Gurgaon Haryana 122050 India
Tel. +91 124 47 64 904

sales@adityadisposed.com
www.adityadisposed.com

The smart choice for sharp blades