



## Skin Graft Blade – 158 mm

Compatible with Watson, Cobbett & Braithwaite Skin Graft Handles

### Product Overview

The KIATO Skin Graft Blade (158 mm) is a sterile, single-use medical device designed for the harvesting of split-thickness skin grafts when used with compatible skin graft handles. Manufactured from medical-grade stainless steel and supplied sterile, the blade is produced under ISO 13485 certified quality management systems to support consistent and reliable clinical use.

### Compliance & Use

This device is intended for use by trained healthcare professionals only. The blade is supplied sterile for single use and should not be reused or reprocessed. Certification details and documentation are available upon request for distributor qualification and tenders.

### Product Specifications

Product Name	Skin Graft Blade
Brand	KIATO
Product Reference	12400
Intended Use	Harvesting of split-thickness skin grafts
Usage Type	Single Use / Disposable
Blade Material	Medical-grade Stainless Steel
Blade Length	158 mm
Blade Thickness	Nominal 0.25 mm
Blade Edge Design	Precision-ground, symmetric double-facet cutting edge
Handle Compatibility	Watson, Cobbett and Braithwaite Skin Graft Handles
Sterility	Supplied Sterile (Gamma Sterilised)
Manufacturing Standard	Manufactured under ISO 13485 certified quality management system
Packaging	Individually sterile packed in easy-peel packs
Shelf Life	5 years from date of manufacture
Export Markets	Global

### Clinical Applications

Used by trained healthcare professionals for split-thickness skin graft harvesting in burn wound management, reconstructive and plastic surgery, trauma-related skin loss, chronic wound management, and post-oncological skin reconstruction. **Procedure selection and technique are determined by the clinician.**



## Product Image & Dimensions



*Illustrative representation for reference only.*

Dimension	Specification
Blade Length	158 mm
Blade Thickness	Nominal 0.25 mm
Handle Compatibility	Watson, Cobbett & Braithwaite



MDSS GmbH  
Schiffgraben 41  
30175 Hannover  
Germany



MDSS CH GmbH  
Laurenzenvorstadt 61  
5000 Aarau  
Switzerland



Aditya Dispomed Products Pvt. Ltd.  
Plot No. 19 Sector 6 IMT Manesar  
Gurgaon Haryana I22050 India  
Tel. +91 124 47 64 904

sales@adityadispmomed.com  
www.adityadispmomed.com

The smart choice for sharp blades